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The screenshot shows the MultiVu software interface with several components labeled:

- Menu Bar**: Located at the top left, containing File, View, Sample, Sequence, Measure, Graph, System, Utilities, Window, and Help.
- Tool Bar**: Located below the menu bar, containing various icons for file operations and control.
- Control Center**: Located below the tool bar, containing buttons for Run, Pause, Abort, and Lock.
- Sequence Command Bar**: Located at the top right, containing a list of sequence commands such as Beep, Call Sequence or Script, Chamber Operations, Remark, Scan Field, Scan Temperature, Scan Time, Sequence Message, Set Field, Set Temperature, Standby, Wait, and Measurement Commands.
- Sequence Window**: Located in the upper left, showing the selected sequence (Sequence1.seq) and its status (Sequence Idle). It also displays the sequence steps: Selected Line 2, Scan Temp from 1.6K to 400K at 1.2K/min in 2 steps, Uniform, Fast Wait For Temperature, Delay 0 secs, No Action; End Sequence.
- Data Window (graph view)**: Located in the lower left, showing a plot of Resistance (Ohm) vs Temperature (K). The plot shows a red curve with data points and error bars, indicating a transition around 2.5 K.
- Data Window (record view)**: Located in the lower right, showing a table of field names and values. The table is as follows:

	Field Name	Field Value
1	Time Stamp (seconds relch)	5090.361
2	Temperature (K)	3.52983522415161
3	Field (Oe)	39.9323474121094
4	Resistance Ch2 (Ohms)	8.929399E-09
5	Resistance Std. Dev. Ch2 (K)	2.320396E-07
6	Phase Angle Ch2 (deg)	70.64775
7	Frequency Ch2 (Hz)	21.3623
8	Averaging Time Ch2 (s)	4.952011
9	AC Current Ch2 (mA)	100
10	DC Current Ch2 (mA)	0
11	Voltage Ampl. Ch2 (V)	2.6946E-08
12	In Phase Voltage Ampl Ch2 (V)	9.929227E-09
13	Quadrature Voltage Ch2 (V)	2.542352E-08
14	Gain Ch2	90000
15	2nd Harmonic Ch2 (V)	4.029424E-05
16	3rd Harmonic Ch2 (V)	0
17	ETO Status Code	411440512
18	ETO Measurement Mode	0
- Status Bar**: Located at the bottom, showing various system parameters such as Sequence Idle, 246.47 K Stable, 0.11 Oe Hold, 0.00 mTol, Tank: 643 Torr, and others.

Figure 3-2. MultiVu Main Window.

PPMS软件介绍—— Menu

File Edit View Sample Sequence Measure Graph Instrument Utilities Window Help

File Edit View Sample Sequence Measure Graph

- New Sequence Ctrl+N
- Open... Ctrl+O >
- Close
- Save Ctrl+S
- Save As...
- User Configuration...
- Print... Ctrl+P
- Print Preview...
- Print Setup...
- 1 T153MT-20240426.dat
- 2 D:\LinFei Liu\20240423.dat
- 3 T153MT-20240425.dat
- 4 D:\LinFei Liu\...\20240418.dat
- 5 eto_default.dat
- 6 ETO_00001.dat
- 7 D:\LinFei Liu\...\20240422.dat
- 8 D:\Guorui Chen\LZX\ETO-.seq
- Exit

File Edit View Sample Sequ

- Undo Ctrl-Z
- Cut Ctrl-X
- Copy Ctrl-C
- Paste Ctrl-V
- Delete Del
- Enable
- Disable

File Edit View Sample Sequence Measure

- Status Bar >
- Control Center
- Tool Bar
- Sequence Command Bar

File Edit View Sample Sequence M

Install...

File Edit View Sample Sequence Measure Graph

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






The screenshot displays the PPMS software interface. It features several control panels:

- Temperature - System:** Shows current Temp (300.00 K) and State (Stable). Control options include Set Point (300.00 K), Rate (5.00 K/min), and Mode (Fast Settle).
- Field:** Shows current Field (-0.30 Oe) and State (Holding). Control options include Set Point (0.00 Oe), Rate (100.00 Oe/sec), and Approach (Oscillate).
- Chamber:** Shows current Pressure (10.0 Torr) and State (Purged). Control buttons include Seal, Purge/Seal, Vent/Seal, Pump Cont., Vent Cont., and HVac.
- Cryostat:** A table showing various parameters: Cryocooler (Running), Magnet (4.32 K, Holding), 4K Plate (4.22 K, 0.84 W), Shield (42.48 K), Loop (673.95 Torr, Circulating, Liquid Full), and Case (Purged).

 A vertical toolbar on the right contains buttons for Refresh Liquid..., System Standby..., Warm Up..., Shutdown..., Restart Cryostat..., and Rapid Shutdown... The status bar at the bottom provides a summary of the system's current state for each component.

Sequence:

- status
- file name
- execution time

Temperature:

- sample temp and status
- setpoint
- set rate and mode

Magnetic Field:

- field and status
- setpoint
- set rate and mode

Sample Chamber:

- pressure
- state

Cryostat and Cryocooler:

- circulation loop pressure
- status
- magnet temp.

Sequence Commands:

System Commands

- Beep
- Call Sequence or Script
- Chamber Operations
- Remark
- Scan Field
- Scan Temperature
- Scan Time
- Sequence Message
- Set Field
- Set Temperature
- Standby
- Wait

Measurement Commands

- Log Data
- Sigma Log Data

Beep – Causes the computer to make the Windows default beep sound.

Call Sequence – Suspends execution of parent sequence file and begins execution of selected child sequence file or script file. When execution of child sequence file or script is complete, execution of parent sequence will continue with the next line.

Chamber Operations – Changes the state of the sample chamber atmosphere (Figure 3-1). See Section 4.5.1 for a description of chamber operations.

Remark – Serves as a message, comment, or visual break for the user only. Does nothing during sequence execution.

Scan Field – Creates a program loop for executing repeated commands at user-defined magnetic field increments. All commands between the “Scan Field...” line in the sequence and the “End Scan” line in the sequence will be repeated at each magnetic field specified by the scan field command. Set the initial and final fields and the scale on which the field increments should appear uniform (linear, H^2 , $H^{1/2}$, $1/H$, $\log(H)$). Also set the total number of field steps. (For uniform linear spacing, you may alternatively set the field increment.) Finally, specify the rate and approach mode used by the magnet controller to achieve each set point:

Scan Temperature – Creates a program loop for executing repeated commands at user-defined temperature increments. All commands between the “Scan Temperature...” line in the sequence and the “End Scan” line in the sequence will be repeated at each temperature specified by the scan temperature command. Set the initial and final temperatures and the scale on which the spacing of the temperature steps should appear uniform (linear, $1/T$, $\log(T)$.) Also set the total number of temperature steps. (For uniform linear spacing, you may alternatively set the temperature increment.) Finally, specify the rate and approach mode used by the temperature controller to achieve each set point:

Sequence Commands:

System Commands

- Beep
- Call Sequence or Script
- Chamber Operations
- Remark
- Scan Field
- Scan Temperature
- Scan Time
- Sequence Message
- Set Field
- Set Temperature
- Standby
- Wait

Measurement Commands

- Log Data
- Sigma Log Data

Scan Time – Creates a program loop for executing repeated commands at user-defined time increments (or immediate repetitions with no time increment.) All commands between the “Scan Time...” line in the sequence and the “End Scan” line in the sequence will be repeated at each time specified by the scan time command. Set the total time in seconds and specify whether the spacing of events should be uniform in time or logarithmic in time. And specify the number of steps. This is the number of times the loop will be repeated. (For uniform spacing in time, you may alternatively specify the time increment.) If the total time is set to zero seconds, then the number of steps defines how many times the loop will be repeated in rapid succession.

Sequence Message – Displays a message on the computer screen and pauses sequence execution until the message is acknowledged or until a timer expires. If the computer is set up with network access and access to a mail server, a message can also be emailed with attachments such as data files.

Set Field – Sets the instrument’s magnetic field. Specify the field, the charging rate, and the approach:

Set Temperature – Sets the sample temperature. Specify the temperature, the rate, and the mode:


Shutdown – Places the instrument in shutdown mode to conserve resources. The field is set to zero, the magnet controller is turned off, the chamber heaters are turned off, and the chamber is sealed.

Wait – Waits for specified conditions to be achieved, then delays a specified amount of time before continuing with sequence execution. Conditions that can be specified to wait for are temperature stability, field stability, and chamber state. This command is usually used immediately after another command in order to make sure the desired outcome of the first command is achieved before proceeding. For example:

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Select one X-axis and up to four Y-axes

Specify the data to plot on each axis

Click to apply additional labeling and scaling to axes

Check for logarithmic axis scaling

Select auto scaling for each axis or specify the min. and max. value for each axis

Axis	Log	Data	Auto	Min	Max
X	<input type="checkbox"/>	Magnetic Field (Oe)	<input checked="" type="checkbox"/>	-100000	250000
Y1	<input type="checkbox"/>	Moment (emu)	<input type="checkbox"/>	-0.6	0.8
Y2	<input checked="" type="checkbox"/>	M. Std. Err. (emu)	<input checked="" type="checkbox"/>	1E-010	0.0001
	<input type="checkbox"/>		<input checked="" type="checkbox"/>	0	0
	<input type="checkbox"/>		<input checked="" type="checkbox"/>	0	0

NOTE: Log scales show only positive values

Data Group: All

Show Time In: Minutes Seconds

Time Display: Relative Absolute

Display Range: All Records Last 539 records From # To # 4539 Prev 4539 Next Records

Filter...

OK Apply Cancel

Display time data as absolute time or relative to first data record in data file in minutes or seconds

Limit the data records plotted to the last few records (like a chart recorder) or to a range of specified records

Click to limit plotted data based on certain data values

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